

The first letter for buyers for BT2-Z:

We sincerely appreciate your purchase of the CoCo Smartwatch BT2-Z. Attached, please find the CoCo BT2-Z User Manual. Please spend a few minutes watching the "Get Started" video listed in the user manual. If you encounter any further issues or need assistance, please refer to the contact information provided within.

Please follow these initial steps to get started with your new BT2-Z:

1. Charge Your Watch: Connect the charging cable to the smartwatch and plug it into a power source. Allow the watch to charge fully before first use.
2. Download the Free CoCo App: Download the CoCo app from the Apple App Store or Google Play Store by searching "CoCo Semcorel."

The app

will help you set up and manage your smartwatch.

3. Register a Free CoCo ID: Open the CoCo app, click "Sign up" at the bottom, and follow the instructions. Ensure your phone number and email are correct; they are used for credentials, health data, and emergency contacts, not marketing.
4. Pair Your Device: Open the CoCo app and follow the first-time user steps to set up a basic profile, add your emergency contacts, and pair your BT2-Z with your smartphone. This will sync your language and time format from your smartphone to your device.
5. Explore Features: Once paired, take some time to explore the various features of your BT2-Z, such as health monitoring, Emergency Care Team (ECT), Medication Reminder, and customizable watch faces.